BRAUN

ThermoScan® **4+ connect**

Ear thermometer

IRT3575

English	4
Spanish.....	24

Contact Us

Questions or Comments?

Call us toll-free at: 1-800-327-7226

Or visit our website at: www.BraunHealthcare.com

Please specify the model number when contacting consumer relations.

Contáctenos

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Llámenos lada gratuita al: 1-800-327-7226

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Kaz USA, Inc., una Empresa de Helen of Troy

400 Donald Lynch Blvd., Suite 300, Marlborough, MA 01752

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ThermoScan® and Braun Family Care™ are trademarks of Helen of Troy Limited.

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

The Bluetooth® word mark and logos are registered trademarks of Bluetooth SIG, Inc. and used under license by Kaz USA, Inc.

App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google LLC.


ThermoScan® y Braun Family Care™ son marcas registradas de Helen of Troy Limited.

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Indications for Use / Intended Use

The Braun ThermoScan® 4+ connect Ear thermometer is intended for the intermittent measurement of human body temperature in a home-use environment.

Warnings and precautions

Keep out of reach of children under 12 years.

NEVER use the thermometer for purposes other than those for which it has been intended. Please follow the general safety precautions when using on children.

NEVER immerse the thermometer into water or other liquids (not waterproof). For cleaning and disinfecting please follow the instructions in the "Care and cleaning" section.

The operating ambient temperature range for this thermometer is 50 – 104 °F (10 – 40 °C).

DO NOT store this thermometer in temperature below -13 °F or over 131 °F (below -25 °C or over 55 °C) or in excessive humidity (above 95% non-condensing relative humidity).

If the thermometer is stored in a location that is cooler or warmer than the room where it will be used, allow it to acclimate to the room temperature for 10 minutes before taking a measurement.

DO NOT use the thermometer if there are signs of damage on the probe or on the thermometer itself. If damaged, **DO NOT** attempt to repair the product.

NEVER insert a sharp object into the probe area or any other open surface on the thermometer.

This thermometer consists of high-quality precision parts. **DO NOT** drop the instrument. Protect it from severe impact and shock.

This thermometer must only be used with genuine Braun ThermoScan® disposable lens filters (LF40/LF20). Other lens filters can affect accuracy.

DO NOT modify the equipment without the manufacturer's authorization.

To ensure accurate measurement **ALWAYS** use this thermometer with a new, clean, disposable lens filter before every measurement.

If the thermometer is accidentally used without a disposable lens filter, to prevent cross contamination, clean the lens after use (see "Care and cleaning").

Keep disposable lens filters out of reach of children.

This thermometer is intended for household use only.

Use of this thermometer is not intended as a substitute for consultation with your physician.

Temperature elevation may signal a serious illness, especially in neonates and infants, or in adults who are old, frail, or have a weakened immune system. Please seek professional advice immediately when a temperature elevation occurs on persons that are:

- Neonates and infants under 3 months (Consult your physician immediately if the temperature exceeds 99.4 °F [37.4 °C])
- Individuals having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer chemotherapy, chronic immunosuppressant treatment, splenectomy)
- Individuals who are bedridden (e.g., nursing home patient, stroke, chronic illness, paraplegia, quadriplegia, surgical recovery)
- A transplant recipient (e.g., liver, heart, lung, kidney)

This thermometer is not intended for pre-term babies or small-for-gestational age newborn babies (birth to 4 weeks of age). Pre-term is defined as babies born before 37 weeks of gestation. Small-for-gestational-age is defined as a newborn baby, born at 37 weeks or later, with a weight below the 10th percentile for newborn babies of the same gestational age.

This thermometer is not intended to interpret hypothermic temperatures. If the device displays a temperature of 97.5°F (36.4°C) or less with an LED backlight color of green or white, and the individual is exhibiting atypical symptoms or behaviors, contact your physician or health care professional.

DO NOT leave children unattended when obtaining temperature measurement.

Please consult your physician if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.

Please consult your physician if the thermometer shows elevated temperature (yellow or red background lights).

Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.

This thermometer contains small parts that can be swallowed or produce a choking hazard to children.

Contains coin battery:

- **INGESTION HAZARD:** This product contains a button cell or coin battery.
- **DEATH** or serious injury can occur if ingested.
- A swallowed button cell or coin battery can cause Internal Chemical Burns in as little as 2 hours.
- **KEEP** new and used batteries **OUT OF REACH of CHILDREN**
- Seek immediate medical attention if a battery is suspected to be swallowed or inserted inside any part of the body.


Pairing to your smart device

Scan the QR code and download the Braun Family Care™ app.


Follow the in-app instructions to pair your ThermoScan® 4+ connect to your smart device.


1. Make sure the device is turned **OFF**. Press and hold the Power button until the screen turns yellow and the Bluetooth icon and the letters "bt" flash. Follow the instructions in the app to finish pairing.


2. When the thermometer successfully pairs to your smart device, the screen will turn green and the Bluetooth icon and letters "bt" will stop flashing. After 5 seconds, the thermometer will be ready to take a temperature.

3. If there is a problem pairing to your smart device, the screen will turn red and Err will be displayed. Follow the instructions in the app to try again.

When the thermometer is actively connected to your smart device the bluetooth icon will show on screen.

Color-coded fever guidance

Color-coded fever guidance helps you to better understand the meaning of your child's temperature with the color indicated on the display. The color-coded screen displays GREEN for a normal temperature, YELLOW for a fever, and RED for a high fever.

Color range	Reading	Meaning
Green	>96.3 °F - 99.4 °F (>35.7 °C - 37.4 °C)	Normal
Yellow	>99.4 °F - 101.3 °F (>37.4 °C - 38.5 °C)	Fever
Red	>101.3 °F (>38.5 °C)	High fever

The Braun Family Care app is able to provide more accurate fever guidance based on age. **The color of the fever level displayed in the Braun Family Care app may not match the color on your thermometer display.** See the FAQ section in the app for more information.

Age range of subject	Green Normal temperature $\geq$	Yellow Elevated temperature	Red High temperature
 0*-3 Months	$\geq 96.4 - \leq 99.4$ °F ($\geq 35.8 - \leq 37.4$ °C)		>99.4 °F (>37.4 °C)
 3-36 Months	$\geq 95.7 - \leq 99.6$ °F ($\geq 35.4 - \leq 37.6$ °C)	>99.6 - ≤ 101.3 °F (>37.6 - ≤ 38.5 °C)	>101.3 °F (>38.5 °C)
 36 Months-adult	$\geq 95.7 - \leq 99.9$ °F ($\geq 35.4 - \leq 37.7$ °C)	>99.9 - ≤ 103.0 °F (>37.7 - ≤ 39.4 °C)	>103.0 °F (>39.4 °C)

*0 months is defined as having a gestational age of greater than or equal to 37 calendar weeks at birth.

Body temperature

Normal body temperature is a range. It varies by site of measurement, and it tends to decrease with age. It also varies from person to person and fluctuates throughout the day. Therefore, it is important to determine normal temperature ranges.

How to use your Braun ThermoScan®

Before taking the temperature, make sure the ear is free from obstructions or excess earwax.

DO NOT attempt to take the reading without a hygienic lens filter on – the measuring sensor must be kept protected and clean from earwax to give accurate readings.


Change the hygienic lens filter with each use to ensure hygiene, to avoid cross-contamination and to ensure earwax residues on the hygienic lens filter do not impact accuracy of the reading.

1. To turn the thermometer on, push the “start” button **1** as shown on page 4-5.

During an internal self-check, the display shows all segments.

The reading shown on the screen is the last temperature taken.

Flashing symbol of lens filter will appear if no lens filter is attached. Temperature can be taken only if lens filter is mounted on the probe.


The thermometer is ready to use after two beeps and the non-flashing ear icon appears.


2. Fit the probe snugly into the ear canal, then push the “start” button.

One beep signals normal temperature. Audio temperature indicator will trigger 2 consecutive beeps when the reading is equal or higher to 99.5 °F (37.5 °C) and 4 beeps if the reading is equal or higher to 101.4 °F (38.6 °C). The result is shown on the display.


3. Gently pull off the hygienic lens filter.

DO NOT use force in removing the hygienic lens filter.

- Put on a new, clean, hygienic lens filter by snapping into position.
- Place the protective cap on.

The Braun ThermoScan® 4+ connect Ear thermometer turns off automatically after 60 seconds of inactivity. The thermometer can be turned on manually by pressing and holding the start button for more than 5 seconds.

Hygienic lens filters (LF40/LF20)

Use a clean hygienic lens filter with each measurement:

- To ensure hygienic use
- To protect from cross-contamination
- To ensure earwax residues do not impact accuracy of the reading

Earwax and dirt on the tip of the hygienic lens filter, even a microscopic layer, can create an obstacle between the sensor and the infrared heat emitted by the ear and impact the accuracy of the measurement.

Temperature taking tips

- Replace the hygienic lens filter with every use to maintain hygiene and accuracy.
- There can be a normal disparity in temperature measurements obtained from the left and right ear. To help ensure that any differences in temperature measurements reflect true changes in body temperature, always take temperature measurements in the same ear.
- The ear must be free from obstructions or excess earwax build-up to take an accurate reading.

- External factors may influence ear temperatures, including when an individual has:
 - been lying on one ear or the other
 - had their ears covered
 - been exposed to very hot or very cold temperatures
 - been recently swimming or bathing, or
 - been wearing hearing aids or headphones for an extended period of time

In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.

- If ear drops or other ear medications have been placed in the ear canal of one ear, take the temperature in the untreated ear.
- Nursing mothers should not take a temperature measurement on themselves or have their temperature measured by another individual while or immediately after nursing a baby, as nursing can affect body temperature.


Care and cleaning

The probe tip is the most delicate part of the thermometer. It has to be clean and intact to ensure accurate readings.

If the thermometer is accidentally used without a hygienic lens filter, clean the probe tip as follows:

Very gently wipe the surface with a cotton swab or soft cloth moistened with an alcohol solution of no more than 70% concentration.

DO NOT use cleaners other than alcohol on the probe tip. Use of other cleaners could damage the tip.

Once the alcohol has dried completely, a new hygienic lens filter can be put on and a temperature measurement taken. If the probe tip is damaged, please contact Consumer Relations.

Use a soft, dry cloth to clean the thermometer display and product exterior. **DO NOT** use abrasive cleaners. **NEVER** submerge this thermometer in water or any other liquid.

Store the thermometer and hygienic lens filters in a dry location free from dust and contamination and away from direct sunlight.




Additional hygienic lens filters (LF40 /LF20) are available at most stores selling the Braun ThermoScan® thermometers.

Troubleshooting

Situation	Solution
 No hygienic lens filter is attached.	Attach new, clean, hygienic lens filter.
 Measurement before device stabilization.	Wait until all the icons stop flashing.
 Ambient temperature is not within the allowed operating range (50–104 °F or 10–40 °C).	Allow the thermometer to remain for 30 minutes in a room where the temperature is between 50 and 104 °F or 10 and 40 °C.
  HI = too high LO = too low Temperature taken is not within typical human temperature range (93.2–108 °F or 34–42.2 °C).	Make sure the probe tip and lens are clean and a new, clean, disposable lens filter is attached. Make sure the thermometer is properly inserted. Then, take a new temperature.
 or  System error.	Reset the thermometer by removing the battery and then replacing it again, wait 1 minute, then turn on again. Please contact Consumer Relations (See Guarantee Card).
 Battery is low.	Insert new battery.
 or  Battery is too low. Device cannot be powered on to the ready state.	Insert new battery.
 Unsuccessful connection.	Check issue in the app for more details.

Changing the temperature scale


1. Ensure the thermometer is **turned ON** and the **probe cover is OFF/REMOVED**.
2. Press and hold down the “start” button for a complete 5 seconds.
3. Release the start button when the desired scale (°F or °C) is shown. This will reset the thermometer to the new temperature scale.

Replacing the battery

The thermometer is supplied with one 3V Lithium cell battery (type CR2032). Insert new battery when the low battery symbol appears on the display.

The battery door design is intended to prevent a small child from easily opening the battery door to address a possible choke hazard. Open the battery compartment by inserting a narrow pointed object such as a ball point pen or paperclip into the small hole above the battery door. While pressing into the hole slowly slide back the battery door with the other hand.

To close battery door align the tabs along the side of the battery into the battery compartment and slowly slide back into position until the door clicks into its locked position.

Only discard empty batteries. To protect the environment, dispose of empty batteries at your retail store or at appropriate collection sites according to national or local regulation

Do not mix alkaline, standard, or rechargeable batteries.

Calibration

This product has been calibrated at the time of manufacture. If used according to the instructions for use, periodic calibration is not required. If at any time you question the accuracy of temperature measurements, please contact Consumer Relations.

FCC warning

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:


- A. Reorient or relocate the receiving antenna.
- B. Increase the separation between the equipment and receiver.
- C. Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- D. Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance can void the user's authority to operate the equipment.

Cybersecurity Guidelines

1. No Anti-Virus or firewall software is used in this device. Keep operating software on your phone up to date.
2. Any internet based communications are authenticated, encrypted, and verified for integrity.
3. Complete backup and restoration of features and procedures to regain configurations are not needed. Critical data points can be restored from the cloud excluding device settings; these will be set via the device set up process.
4. In the event of a device reset, no recovery of device settings is possible.
5. The Mobile App needs Bluetooth and Location permissions for the thermometer to operate as intended with the app.

6. No supplemental secure configuration is required by the user. However, we encourage end-user to enable user-authentication in their smartphone such as pin code, swipe or biometric authentication like face recognition or fingerprint recognition.
7. App Updates when needed will be pushed through notifications to the smartphone mobile device and user will choose when to update. Software updates to the thermometer will be pushed through the App and user must follow the provided instructions to begin the update.
8. In the event of an anomalous connection to the network or loss of connection to the network, the device will continue to operate as intended.
9. Settings for the thermometer are not retained in the cloud or recovery.

11. For software components please contact Customer Service. Contact information is on page 2.
12. There are no user configurable connections. If you have questions regarding cyber-security questions, please contact customer support at the number provided.
13. Helen of Troy expects to maintain all online services and Apps for the ongoing future to ensure your satisfaction. However, in the event that an "end of support" situation becomes unavoidable, followed by a grace period of 3 years and then an "end of life", notification will be provided via the App. In that situation, Helen of Troy will provide a final App and firmware update to allow continued offline use. Following the update, no ongoing support for the device or App will be provided.

3 Year Limited warranty

You should first read all Instructions before attempting to use this product.

- A. This 3 year limited warranty applies to repair or replacement of product found to be defective in material or workmanship. This warranty does not apply to damage resulting from commercial, abusive, unreasonable use or supplemental damage. Defects that are the result of normal wear and tear will not be considered manufacturing defects under this warranty. **KAZ IS NOT LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY NATURE, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THIS PRODUCT IS LIMITED IN DURATION TO THE DURATION OF THIS WARRANTY.** Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you also may have other rights which vary from jurisdiction to jurisdiction. This warranty applies only to the original purchaser of this product from the original date of purchase.
- B. At its option, Kaz will repair or replace this product if it is found to be defective in material or workmanship.
- C. This warranty does not cover damage resulting from any unauthorized attempts to repair or from any use not in accordance with the instruction manual.

Call us toll-free at 1-800-327-7226 or

Email: ConsumerRelations@HelenofTroy.com

Please be sure to specify a model number.

NOTE: IF YOU EXPERIENCE A PROBLEM, PLEASE CONTACT CONSUMER RELATIONS FIRST OR SEE YOUR WARRANTY. DO NOT RETURN THE PRODUCT TO THE ORIGINAL PLACE OF PURCHASE. DO NOT ATTEMPT TO OPEN THE DEVICE HOUSING YOURSELF, DOING SO MAY VOID YOUR WARRANTY AND CAUSE DAMAGE TO THE PRODUCT OR PERSONAL INJURY.

Model:	IRT 3575
Displayed temperature range:	93.2 °F – 108 °F (34 °C – 42.2 °C)
Operating ambient temperature range:	50 °F – 104 °F (10 °C – 40 °C)
Storage temperature:	-25°C - 55°C (-13°F - 131°F)
Display resolution:	0.1°F or 0.1°C
Accuracy for displayed temperature range:	±0.4 °F (±0.2 °C) in range 35.5-42°C (95.9-107.6°F) ±0.5 °F (±0.3 °C) outside that temperature range

Clinical accuracy characteristics and procedures are available upon request.

Measuring Site:	Ear
Reference Site:	Oral (estimated oral temperature is displayed)

Operating and storage

Relative humidity:	15% to 95% non-condensing
Battery life:	3 years /800 measurements
Service life:	5 years
Warranty:	3 years

If device is not used or stored within specified temperature and humidity ranges the technical accuracy cannot be ensured.

This infrared thermometer meets requirements established in ASTM Standard E 1965-98 (for the thermometer system [thermometer with lens filter]). Full responsibility for the conformance of the product to the standard is assumed by Kaz USA, Inc.

ASTM laboratory accuracy requirements in the display range of 37 to 39°C (98 to 102°F) for IR thermometers is ±0.2°C (±0.4°F), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is ±0.1°C (±0.2°F).

This thermometer is specified to operate at a pressure of 1 atmosphere or at altitudes with a pressure of up to 1 atmosphere.

Type BF Applied
Part

See
instructions
for use

Operating
temperature

Hygienic lens filter:
Do not reuse

Storage
temperature

Keep it dry

Manufacturer

Date of
Manufacture

IP22

Protection
provided by
enclosures

Attention, consult
accompanying
documents

Storage
relative
humidity

Operating
relative
humidity

Recycle:
Electronic
Equipment

Subject to change without notice.

This device conforms to the following standards:

ASTM E1965-98 - Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature.

ISO 80601-2-56 – Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement.

ANSI AAMI IEC 60601-1 - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.

IEC 60601-1-2 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances – Requirements and tests.

IEC 60529 - Degrees of protection provided by enclosures (IP Code)

NOTE: DO NOT use this device in the presence of electromagnetic or other interference outside the normal range specified in IEC 60601-1-2.

IEC 60601-1-11 - Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

ISO 10993-1 - Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.

ISO 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements.

Internally Powered Equipment
Continuous Operation
Not Protected against Ingress of Water.

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC. For detailed description of EMC requirements please contact Consumer Relations.

Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.

This product contains batteries and recyclable electronic waste. To protect the environment, do not dispose of it in the household waste, but take it to appropriate local collection points.

Guidance and manufacturer's declaration – electromagnetic immunity

The IRT3575 is intended for use in the electromagnetic environment specified below.
The customer or the user of the IRT3575 should ensure that it is used in such an environment.

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
		HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields ^{a)}	IEC 61000-4-3	10 V/m ^{f)} 80 MHz – 2,7 GHz ^{b)} 80 % AM at 1 kHz ^{c)}
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See RF wireless communication equipment immunity table below
RATED power frequency magnetic fields ^{d) e)}	IEC 61000-4-8	30 A/m ^{g)} 50 Hz or 60 Hz

a) The interface between the PATIENT physiological signal simulation, if used, and the ME EQUIPMENT or ME SYSTEM shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the ME EQUIPMENT or ME SYSTEM.

b) ME EQUIPMENT and ME SYSTEMS that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test.

c) Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

d) Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry.

e) During the test, the ME EQUIPMENT or ME SYSTEM may be powered at any NOMINAL input voltage, but with the same frequency as the test signal (see Table 1).

f) Before modulation is applied.

g) This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that the ME EQUIPMENT or ME SYSTEM will be used closer than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.

Guidance and manufacturer's declaration – electromagnetic emissions

The IRT3575 equipment is intended for use in the electromagnetic environment specified below.
The customer or the user of the IRT3575 should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The ME equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	Complies
Harmonic emissions IEC 61000-3-2	Not Applicable	The ME equipment is solely battery powered.
Voltage fluctuations/ flicker emissions	Not Applicable	


Guidance and manufacturer's declaration – RF wireless communication equipment immunity

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5 500						
5 785						

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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© 2024, All Rights Reserved.